

BATCH RECORD · LOT-2026-100

Product: mAb-X Date: 2026-03-10

Step 1: Buffer Preparation

pH: 7.2

Temperature: 25.0 °C

Volume: 500 mL

Start Time: 08:30 Operator: JKL

Notes: Solution was clear.

Step 2: Centrifugation

Speed: 3000 g

Duration: 15 min

Temperature: 4.0 °C

Start Time: 09:15 End Time: 09:30 Operator: JKL

Notes: Pellet compact, supernatant clear.

Step 3: Filtration

Filter Size: 0.22 μ m

Pressure: 5 PSI

Start Time: 09:45 Operator: JKL

Notes: Filtration completed without issues.