

Multi-Role Buffer Protocol — Batch Record

Run: Preview | Status:

Started: | Completed:

Project: | Organization:

SOP Reference: v2 (effective)

Procedure Execution

Media Prep

Step	Description	Value	Operator	Reviewer	Scheduled	Actual Start	Actual End	Notes
Weigh Reagents	Weigh NaCl and KCl. NaCl: 8, and KCl: 0.2.	NaCl (g): KCl (g):						
Dissolve	Add to 1000 mL water and stir.							

QC

Step	Description	Value	Operator	Reviewer	Scheduled	Actual Start	Actual End	Notes
Measure pH	Measure and adjust pH. Target pH: 7.4.							

Approval