

Date: \_\_\_\_\_

Lot/Batch #: \_\_\_\_\_

## Media Prep

| # | Step           | Description                                                              | Value / Result                 | Initials |
|---|----------------|--------------------------------------------------------------------------|--------------------------------|----------|
| 1 | Weigh Reagents | Weigh NaCl and KCl. NaCl: 8, and KCl: 0.2.<br><i>Notes:</i> See Figure 1 | NaCl (g): 8.01<br>KCl (g): 0.2 | D.S.C.   |
| 2 | Dissolve       | Add to 1000 mL water and stir.<br><i>Notes:</i> See Figure 2             | 1000                           | D.S.C.   |

## QC

| # | Step       | Description                                                          | Value / Result | Initials |
|---|------------|----------------------------------------------------------------------|----------------|----------|
| 1 | Measure pH | Measure and adjust pH. Target pH: 7.4.<br><i>Notes:</i> See Figure 3 | 7.38           | J.W.     |

## **Notes & Observations**

**Dr. Sarah Chen | 2026-04-01 10:35:00 | ACTIVE**

Observed slight turbidity after buffer prep.

**Dr. Sarah Chen | 2026-04-01 11:20:00 | ACTIVE**

[ANOMALY] ANOMALY: Cell viability below threshold at 85%.

## Figures

*Figure 1: test\_figure\_1.png*

Figure 2: test\_figure\_2.png

*Figure 3: test\_figure\_3.png*

## Role Sign-Off

Media Prep

---

Signature / Date

QC

---

Signature / Date