

Multi-Role Buffer Protocol

Protocol with Media Prep and QC roles.

Media Prep

1. Weigh Reagents

Weigh NaCl and KCl.

NaCl: 8, and KCl: 0.2.

Allow 10 minutes for this step.

2. Dissolve

Add to 1000 mL water and stir.

Allow 15 minutes for this step.

QC

1. Measure pH

Measure and adjust pH.

Target pH: 7.4.

Allow 5 minutes for this step.

Approvals

Prepared by (Name / Date)

Reviewed by (Name / Date)