

# BATCH RECORD · LOT-2026-201

Product: CHO-K1 Expansion    Date: 2026-03-12

## Step 1: Cell Seeding

Cell Density: 0.5e6 cells/mL

Media Volume: 50 mL

Passage Number: 12

Operator: ABC

## Step 2: Incubation

Temperature: 37.0 °C

CO2: 5 %

Duration: 72 hours

Operator: ABC

### Step 3: Harvest

Viability: 95 %

Cell Count: 2.1e6 cells/mL

Operator: ABC