

Date: _____

Lot/Batch #: _____

Media Prep

#	Step	Description	Value / Result	Initials
1	Weigh Reagents	Weigh NaCl and KCl. NaCl: 8, and KCl: 0.2. <i>Notes:</i> See Figure 1	NaCl (g): 8.01 KCl (g): 0.2	D.S.C.
2	Dissolve	Add to 1000 mL water and stir. <i>Notes:</i> See Figure 2	1000	D.S.C.

QC

#	Step	Description	Value / Result	Initials
1	Measure pH	Measure and adjust pH. Target pH: 7.4. <i>Notes:</i> See Figure 3	7.38	J.W.

Notes & Observations

Dr. Sarah Chen | 2026-04-01 10:35:00 | ACTIVE

Observed slight turbidity after buffer prep.

Dr. Sarah Chen | 2026-04-01 11:20:00 | ACTIVE

[ANOMALY] ANOMALY: Cell viability below threshold at 85%.

Figures

Figure 1: test_figure_1.png

Figure 2: test_figure_2.png

Figure 3: test_figure_3.png

Role Sign-Off

Media Prep

Signature / Date

QC

Signature / Date